

Clinical Experience with the Statera vi Cordless Endodontic Motor in Routine and Complex Endodontic Cases

Prof. Igor Tsesis, Israel

Abstract

This report outlines my experience using the *Statera vi Cordless Endodontic Motor* over several months in a diverse range of endodontic procedures. The evaluation focuses on the device's ergonomics, functionality, and performance, particularly its adaptive motion technology, integrated apex locator, and overall impact on treatment workflow and patient management.

Introduction

Efficient and predictable root canal instrumentation remains a cornerstone of successful endodontic therapy. The integration of intelligent motor control, user-centric design, and digital enhancements has led to significant advancements in endodontic motors. The *Statera vi Cordless Endodontic Motor* is designed to address clinician demands for precision, adaptability, and streamlined workflow in both routine and anatomically complex cases.

Clinical Experience

Over the course of several months, the *Statera vi* was incorporated into daily clinical use, encompassing a wide range of treatments, from straightforward single-canal cases to more demanding multi-rooted teeth with curved or narrow canals. My impressions are based on handling, motor responsiveness, safety features, and integration into the endodontic protocol.

Ergonomics and Handling

The *Statera vi* has a well-balanced, ergonomic feel. Its lightweight, cordless form factor significantly reduced hand fatigue, particularly during longer procedures. The cordless design improved access, especially in posterior regions and in patients with limited mouth opening, enhancing both clinician comfort and procedural efficiency.

Motor Performance

The motor exhibited quiet and smooth operation with consistent torque delivery. The torque response was predictable, helps maintain control during canal shaping and enhancing both safety and precision.

Balance Force Adaptive Motion Technology

A key feature of the *Statera vi* is its Adaptive Balanced Force capability, which responds dynamically to canal resistance. This was particularly effective in reducing stress on the file system and appeared to mitigate the risk of instrument separation, especially in curved or constricted canals.

Integrated Apex Locator

The built-in apex locator has worked well in my hands, eliminating the need for a separate device and enabling continuous workflow during working length determination. This integration reduced procedural time and equipment handling.

User Interface and Workflow Integration

The device's interface is intuitive and easy to navigate. Switching between file systems or modifying settings was straightforward and did not disrupt clinical flow. The system accommodates various file protocols, making it adaptable to different clinician preferences.

Battery and Charging

Battery performance was more than adequate for daily clinical use. The compact charging base was unobtrusive and integrated well into the operatory setup.

Conclusion

The *Statera vi Cordless Endodontic Motor* has proven to be a highly functional, intuitive, and ergonomic device that supports safe and efficient endodontic therapy. Its balance force adaptive motion technology and integrated apex locator contribute meaningfully to clinical outcomes and workflow efficiency. While long-term durability and performance will become clearer over extended use, the motor has thus far demonstrated strong reliability in clinical practice. It is a valuable addition to the contemporary endodontic armamentarium.

Prof. Igor Tsisis, Israel

February 22, 2025

