




BINGO PRO APEX LOCATOR

- Fully automatic
- No more guess work, ever
- Interactive operation, easy-to-use
- Built-in Demo Mode
- User-friendly intelligent interface
- VIRTUAL APEX function
- Audio feedback with volume control
- Built-in Self-Test and protection
- Auto-off
- Rechargeable battery


BINGO3DPRO APEX LOCATOR



File tracking



Large Apical zone



Over-instrumentation

Specifications

Dimensions:	L 130 x W 80 x H 60 mm
Weight:	300 gr.
Display type:	Color Graphic TFT
Window size:	3.5"
Power supply:	2.4V NiMH rechargeable battery
Charger:	Input: 230V / 50-60 Hz or 120V / 50-60Hz
Output DC:	6V DC
Warranty:	24 months

